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EN

WARNING: HEALTH AND SAFETY INFORMATION

Follow the precautions below to avoid damaging your device, accessories, or any connected devices, and to reduce the risk of personal injury, property damage, discomfort or other potential hazards. Keep all safety instructions and other instructions for future reference. For safety reasons we recommend using only original consteed components. All modifications and changes to the consteed components are prohibited. The manufacturer assumes no liability for any resulting damage.

Target audience

Consteed the smart horse sneaker is designed for Hobby riders, semi-professional sport riders, farriers, physiotherapists or vets tracking and seeking data-driven enhancements to the health, well-being and training status of their horses. To discover how our product can specifically benefit you, please visit our website www.consteed.com.

Consteed Set

Handle your consteed Set with care. You may damage the components by puncturing, burning, deforming, crushing, dropping, modifying, attempting to insert foreign objects into or by disassembling them.

  **WARNING:** Battery overheating or injury may be caused when using a damaged smart horse sneaker, protection mat, or other component of the consteed Set.

The smart horse sneaker

Avoid anything that will affect the operation of the smart horse sneaker. This includes cutting, heating, melting or excessively pulling, stretching, twisting, or bending the shoes. Check that the protection mat in your shoes, as well as the shoes itself on the hoofs fit properly. Otherwise, this could result in a damaged product. Also make sure to use the shoes only in a clean and dirt-free condition. As our consteed sensoric and the upper protection mat should be positioned under the hoof of the horse, make also sure that their undersides and/ or irons are clean and dirt-free as well. Otherwise, this lead to incorrect data registration.

 **INFO:** Always avoid using the shoes without a protective mat to prevent sensor damage!

Electronic Unit

Please make sure, that the area which includes the charging-contact pins is clean and dirt-free.

The electronic unit / control unit is permanently installed with the smart horse sneaker via an electronic unit holder. Do not open the holder or the housing of the electronic unit to prevent damage, to prevent damages. To ensure trouble-free operation, you should avoid knocks, falls and external influences such as mud, puddles of water, deep snow or humidity of over 95%. Always ensure that the electronic unit is on the outside of the hoof boots and hooves to prevent damage when the horse move.

Battery and Charging

Our smart horse sneaker are powered by a lithium battery, which is a sensitive component that can cause injury if damaged.   **WARNING:** Do not attempt to remove the battery. Replacement by unqualified professionals may cause damage to your device. The usage of an unqualified or improper battery may constitute a potential risk of fire, explosion, leakage, or other hazards. For battery service, reach out to the

consteed customer service via the website at www.consteed.com or email address service@consteed.com. The battery can be charged by plugging the smart horse sneaker to the provided magnetic y-charger which must be connected to the provided USB-A adapter or an USB-A port which is able to deliver at least 300 mA. Please only use the enclosed charging y-cable. Never pinch or twist the USB-A charging cable and do not force the adapter to a port.

ⓘ INFO: The smart horse sneaker normal charging voltage is $5\text{ V} \pm 0.25\text{ V}$ with a maximum charging current of 150mA per connected horse sneaker. The maximum input voltage of 5.25 V shall not be exceeded under any circumstances. For your own safety make sure to only use certified charging equipment.

⚠ WARNING: Fire, electric shock, injury, or a damaged Consteed smart horse sneaker or other property can be caused by using damaged cables or chargers. The same applies to charging when moisture is present or due to failure to use compatible charging accessories.

ⓘ INFO: Do not charge the smart horseshoe below 32°F, as this may impair charging or damage the battery.

Radio Frequency and Medical Device Interference

Observe the regulations prohibiting the use of radio technology (e.g. cellular or Wi-Fi). A Bluetooth low energy transceiver operating at 2.4 GHz SRD band is featured in the consteed smart horse sneaker. The maximum transmission output is above the standard required for the end device. Our consteed smart horse sneaker complies with regulations governing radio frequency emissions. Nevertheless, such emissions may have a negative effect on the operation of other electronic equipment. For example, while flying, use the consteed smart horse sneaker product only in accordance with the airline's instructions. In the case that you need to disable the Bluetooth functionality of the consteed smart horse sneaker, follow the instructions of the user manual to turn off the boots.

ⓘ INFO: Some components of the consteed set emit electromagnetic fields. Since they could interfere with defibrillators, pacemakers, or other electronic medical devices, we advise you to consult your physician or your medical device manufacturer to ensure that it is safe to simultaneously use your implanted device and this product. For further questions about using our product near implanted medical devices consult your healthcare provider before using the consteed set. The use of the product is at your own risk.

Medical Device

The consteed set is exclusively designed and created for horse health purposes and is not a medical device. Consult your physician, veterinarian, or physiotherapist as soon as there are uncertainties regarding the measurement results. The information we provide serves not as a basis for decision making and should not be relied upon without using other sources of information, like mentioned before.

⚠ WARNING: Stop measuring immediately in the unlikely event of any slipping, twisting or even delocalization of the shoe, if the horse shows signs of restlessness or if abrasions are detectable on the coat.

Storage and Transport

When not in use, fully charge the smart horse boot sneaker and store disconnected from the smartphone or any other electronic device to extend battery life. The sent along packaging bag is there to protect the consteed smart horse sneaker components. Therefore, insert the components free of dirt and dry into the transport bag and check if it is fully closed to withstand environmental influences. In general, we recommend storing the components at room temperature and to not expose the set to direct sunlight.

ⓘ INFO: Do not store the consteed smart horse sneaker less than 10 degrees and no longer than six months without charging. To avoid a loss of capacity, the battery must be recharged after a maximum of 3 months of storage. When not in use, the relative humidity is fixed at $60 \pm 20\%$.

Cleaning and Drying

The product needs to be dry and clean to operate in its best possible way. We recommend using a slightly damp cloth or a brush to do so. Remove the protection mat for cleaning so that you can also remove any dirt underneath. **⚠️ WARNING:** Do not use harsh chemicals or detergents under any circumstances. Cleaning under abundant water, extreme temperature or in the washing machine/dish washer is strictly prohibited as well. After using the consteed the smart horse sneaker in damp conditions, it is very important to remove the protection mat and let all components to dry naturally. A well aired room is most suited for this purpose. Make sure to stay away from non-natural drying/heating equipment or direct sunlight. This applies for all components of the consteed set.

Environmental Characteristics

Operating temperature and humidity:

The consteed set ideally operates in temperatures between -41°F and 113 °F. The max. relative humidity is fixed at 95%. For charging the batteries the min. temperature is 32°F The maximum temperature is 104 °F.

Distraction

Using our product while performing certain activities may distract you. To reduce the risk of accidents when using the smart horse sneaker, especially when handling the horse, use the product on an obstacle-free, non-slip and level surface. In addition, the animal should be in a familiar and quiet environment with plenty of free space to avoid unpredictable disturbances and noises and to prevent the risk of crushing. It must not be used while engaging in other activities like operating heavy machinery, operating vehicles, etc. Additionally, always think of the horse's welfare when using the device. Furthermore, local laws about using mobile devices and headphones must be followed.

Handling horses

The consteed set should only be used on horses by people who are experienced and sensitive to the nature of the animal. To avoid injury to humans and animals, proceed with caution when leading the horse, have a calming effect on the horse by speaking, involve a helper during the measurement if possible and position the horse's stance and line of sight to the outside to calm the horse and demonstrate an escape option.

Water Resistance

The consteed set is splash-proof, but not fully waterproof. Do not let it submerge. Do not wash the product and clean it only according to the instructions provided in this user manual. A damaged product will result in compromised humidity resistance.

Age Limitation

This product is not a toy. It must be kept out of reach of children. The usage by children is not intended.

⚠️ WARNING: Some product components contain small parts that are not suitable for children under the age of three years. If swallowed, a doctor must be contacted immediately. There is also a risk of strangulation for children with the enclosed Y-charging cable. For purchasing and using our product, it is required by law that you confirm you are at least or over 18 years of age or have parental consent.

Repair and Service

If the product is not working properly or has been damaged, you are advised to contact customer service for authorized service. **⚠️ WARNING:** Do not attempt to repair the consteed set yourself. Disassembling the product may result in loss of humidity resistance, accuracy, general functionality and other product properties.

Injury may be caused by disassembling. Unauthorized changes or any alterations to the device will lead to loss of warranty coverage.

EU/UK-Compliance

Konrad Hornschuch AG hereby confirms that this product complies with the fundamental requirements and other relevant provisions of Directive 2014/53/EU and comparable EU directives.

 ^{UK} For further information, please head to the declaration of conformity section at www.consteed.com.

USA-Compliance (coming soon)

This device complies with part 15 of the FCC rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Users are not permitted to make changes or modify the device in any way. Changes or modifications not expressly approved by BAL.ON will void the user's authority to operate the equipment.

This device and its antenna must not be co-located or operating in conjunction with any other antenna or transmitter. Non-body-worn devices must be placed at least 8" (20cm) away from the body.

Waste Electrical and Electronic Equipment (WEEE) Directive and Battery Directive

The consteed smart horse sneaker contain a lithium-ion battery. The WEEE symbol states that only qualified recyclers or hazardous material handlers are allowed to dispose the product when the lifespan reaches its end. This battery must not be mixed with the solid waste stream.

 **WARNING:** Improper disposal may lead to fire, explosion, and/or other hazards. Recycling our product and a separate collection of its electrical accessories and battery will help to conserve natural resources, the environment and protect human health.

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Konrad Hornschuch AG, Salinenstrasse 1, 74679 Weissbach, Germany.
Konrad Hornschuch AG is part of the Continental corporation.